

# **MOTION REDEFINED**



  
SIGNATURE<sup>®</sup> DUO  
Dual Mobility Cup

*Made in UK*

# SIGNATURE<sup>®</sup> DUO


## Signature DUO Dual Mobility Cup

**Signature DUO Dual Mobility System** by Biorad Medisys is designed to provide increased range of movement and to decrease the risk of dislocation. It is a prominent solution to both primary and revision **THA** reducing wear and enhancing stability without sacrificing natural movement


### Benefits

- Reduced risk of dislocation
- Higher range of motion
- Reduced wear rate


### Combining the best of two worlds




## Signature DUO- Acetabular Cup




## Signature DUO- Neutral Liner



## Signature DUO- Acetabular Cup



## Signature DUO- Neutral Liner



## Compatibility chart

Signature DUO Cup Porous coated Cup	Signature DUO XLPE Poly Insert	CoCr Femoral Metal Head	
		22.2	28
44	44/22.2-46/22.2	Yes	—
46			
48	48/28-50/28	—	Yes
50			
52	52/28-54/28	—	Yes
54			
56	56/28-58/28	—	Yes
58			
60	60/28-62/28	—	Yes
62			



BRP/CAT/OR/SDU/01/00: 06-23