



SIGNATURE[®] NX
UNCEMENTED HIP SYSTEM

Fixation Redefined




Trabecular Structure


The next best thing after bone

SIGNATURE NX CUP

The trabecular acetabular cup by **Biorad Medisys** is the **best choice after bone**. It is made from titanium alloy for use in all kinds of hip replacement surgeries (primary and revision hip arthroplasty).



THE TRABECULAR STRUCTURE



The Cup of Choice

for all primary and revision hip replacement surgeries



Signature NX Cluster Hole Cup



Signature NX Multi-Hole Cup

SIGNATURE NX CUP

Material Specs: Ti6Al4V ELI

Liner Cup (mm): 44, 46, 48, 50, 52, 54, 56, 58, 60, 62, 64

BIOLOX® delta HEAD

Material Specs: Fourth Generation Advanced Ceramic
Lowest wear rate
Outstanding biocompatibility
Excellent stability



Sizes (mm): 28, 32, 36

LINER

Material Specs: Highly Cross Linked Polyethylene(XLPE)

Wears less & lasts longer

Option of neutral & 20° hooded liner

Star-shaped profile to avoid rotational movement



Cup Sizes	44	46	48	50	52	54	56	58	60
Neutral/ Hooded Liner	28		28/32		32/36				

Sizes in mm

CLUSTER-HOLE

Sr.no	Diameter (mm)
1	44
2	46
3	48
4	50
5	52
6	54
7	56
8	58
9	60
10	62
11	64

MULTI-HOLE

Sr.no	Cup Size (mm)	No.Of Holes
1	44	6
2	46	8
3	48	8
4	50	8
5	52	10
6	54	10
7	56	10
8	58	10
9	60	10

bioradmedisys™
science for people

Corporate Office
Biorad Medisys Pvt. Ltd., Survey No.: 48/3 & 48/7,
Pashan Sus Rd, Sus Village, Taluka Mulshi, Pune - 411021. MH, India.
P: 020-35000200 | M: 8605422000
bioradmedisys.com

Mumbai Office
Biorad Medisys Pvt. Ltd., 221, Damji Shamji Industrial complex,
9 L.B.S. Marg, Kurla (West), Mumbai-400 070. India
P: +91-22-2502 1474/3704/1546

Follow us on    